



USER GUIDE

WWW.HEALTHYWAVEMAT.COM

Healthy Wave Mats offers the world's widest variety of affordable high-quality far infrared * PEMF * Red light * Negative Ion * Natural crystals mats.





CONTENT

Quick Start Guide	4
Usage and Tips	4
Warnings	5
Temperature Control.	6
LED Controller Instructions	7
Comfort	11
Jet Controller Instructions.	12
Platinum Controller Instructions	14
Care	16
Notes	17
Featured Therapies.	18
Our Policies	20
Troubleshooting	21
Disclaimer.	22
Frequently Asked Questions	23



THANK YOU

At Healthy Wave Mat, we strive to achieve complete customer satisfaction by providing high-quality products that will exceed your expectations. Our products incorporate the most effective natural materials and technologies that reach the common goal of improved health. It's very important for us to make a difference in our customers' lives, which is why we offer up to a 5-year limited warranty for our products.



QUICK START



1

CONTENT REVIEW

Make sure you have your controller and warranty card



2

PLUG IN

Plug one end of the controller into the mat



3

POWER UP

Plug the other end into a surge protector



4

POWER ON

Press the Power button. Refer to the Controller Instructions (see page 7)



USAGE AND TIPS

1. Please follow Safety Precautions at all times when using your heating pad (mat).
2. Preheat your mat for 15-20 min before use.
3. We suggest starting the first few sessions with 104°F (40°C) for 30-minute sessions, a few times a week.
4. Gradually adjust the temperature/time settings until you find a comfortable level. Just pay attention to how you feel at different temperature settings.
5. Only sleep on a mat with low heat settings. Extended use at high temperatures may result in dehydration.
6. You can use your mat every day with low and medium settings. At high temperatures, you should not have more than one session per day.
7. Please allow up to 40 min for the surface to reach the highest temperature.
8. Stay hydrated. If you use the mat at higher temperatures, your body may lose water. It is important to drink plenty of water before and after use.
9. Cover the mat to protect it from moisture and stains, to protect your body from direct contact with hot gemstones, or to alleviate some heat if the surface is too hot for you. Make sure you cover the entire surface. Using a cover also helps gemstones to heat up faster.
10. A cover does not reduce the effects of FIR, negative ions, or PEMF.



Lie on the mat in any way you choose. With low settings you can lie on the bare mat



With medium and high temperature settings you should cover your mat with a towel or an optional cover to protect it from spills, sweat, or stains



WARNINGS

1. Please take time to properly familiarize with the information in this user guide before using your mat. It is the owner's responsibility that all individuals who may use the mat are aware of the correct procedures.
2. Use your mat on flat surfaces. Soft and flexible models can be used on any reasonably flat surface, while firm models should be used on flat surfaces only. A bed, floor, or massage table makes for an ideal spot. Belts, vests, pillows, and other special units are easy to use.
3. When using your mat in a bed or chair, we recommend placing a foil blanket or bath towel between the backside of the mat and the surface of the bed/chair as the backside of the mat may slightly heat up.
4. Protect your mat from moisture at all times. This includes spills and perspiration. You must cover your mat. You may use a provided waterproof cover, towel, or blanket. For a belt, vest, pillow, or other special unit wear clothes that absorb moisture well.
5. When you cover your mat, always cover the entire surface. **Do not partially cover your mat in any way.** Partial coverage can cause the heat sensors to malfunction.
6. Never leave any objects on any part of the mat surface while the mat is heating, including the controller. This can lead to uneven surface heating and possible damage to the mat or heat damage to the object that is placed on the mat.
7. Do not attempt to take apart or open the mat or controller for any reason. If an issue arises, please contact us first.
8. We recommend you always use a surge protector. The wattage of the surge protector should be the same or higher than that of the mat.
9. Never use pins or other metallic fasteners to hold a mat in place.
10. **Do not handle the plug with wet hands.**
11. Keep a controller ventilated to prevent it from overheating. Do not put the controller on soft surfaces that may cover the vents. If the device overheats due to misuse, turn off the heat and place the controller on a table or other firm surface.
12. If you have any preexisting conditions that may be cause for concern, or if you have either heart disease, or hypertension, or any other serious medical condition you should consult with a doctor before using this or any medical device.
13. You should avoid heat therapy if you have vascular disease, deep vein thrombosis, or multiple sclerosis.
14. Limit heat treatment at maximum temperature to one hour per day to prevent dehydration and hyperthermia.
15. Pregnant women, those who are nursing, and individuals who have a pacemaker, should not be exposed to the increased heat from the mat.
16. **PEMF therapy is not suitable for those who are pregnant or individuals with a pacemaker.**
17. Always unplug and power down the mat completely before folding it to put away.
18. Only fold the mat in its predetermined positions.
19. Do not immerse in water
20. **Never unplug the mat by pulling its connecting cord.**
21. Do not use an electrical outlet that has become loose or does not engage the power plug tightly.
22. Do not remove product labels.
23. Do not use the mat with liniment, salve, or ointment preparations that contain heat-producing ingredients. Skin burns can result.
24. Do not leave your mat unattended while it is powered on, especially if children are present.



ATTENTION: Always cover your mat completely.

TEMPERATURE CONTROL

Heating the gemstones in your mat naturally releases far-infrared rays and negative ions. Adjusting the strength of your exposure is as simple as adjusting the temperature of your mat. The actual temperature of a mat may vary depending on the model. The gemstones used are natural heat conductors: they get hot quickly and slowly release their heat into the air. As a general rule of thumb, the higher you set the temperature, the shorter the session should be.

86°F - 95°F (30-35°C) | OVERNIGHT

Good for prolonged and overnight use. Promotes better sleep and the relaxation of both mind and body. The safe setting for those who have a sensitivity to heat.

105°F - 130°F (45-55°C) | UNDER 2 HR


Good for cycles of 2 hours or less*. Temporary relief of joint pain, stiffness, and minor muscle pain; provides muscular relaxation, and a temporary increase in local circulation where applied. Not too warm for the average person.

130°F - 150°F (55-65°C) | 30-90 MIN

Good for 30-90 min cycles*. Temporary relief of joint pain, stiffness, and minor muscle pain; provides muscular relaxation, and a temporary increase in local circulation where applied. Not too warm for the average person. Should only be used twice a day at this setting and nothing more.

150°F - 160°F (65-70°C) | 30-60 MIN

Suitable for 30-60 min cycles*. Provides a far-infrared sauna experience when the heat is trapped. Cover your body with the provided mylar blanket or with a towel to trap the heat. When you relax in gentle far-infrared heat, the body works hard to release sweat, temporarily improving circulation, and recovering faster. Should only be used once a day at this setting and nothing more.



*A mat has to be covered in order to protect your body from direct contact with hot gemstones, and to protect the mat from moisture and stains (see Usage and Tips, #9)

CONTROLLER INSTRUCTIONS

HEAT CONTROLLER


POWER



Press the Power button to turn ON the controller. TEMP and TIMER indicators will light up on the screen.

TEMPERATURE SET & ADJUST

Press the TEMP button once to turn ON the heat. Press TEMP again to initiate temperature adjustment mode (the temperature indicator will start flashing). Set the desired temperature and press TEMP again to save the setting.



°F or °C UNITS SET


Press & Hold
for 6 sec



To switch between
Fahrenheit & Celsius.

TIMER SET & ADJUST

Press the TIME button to initiate time adjustment mode (the time indicator will start flashing). Set the desired time and press the TIME button again to save the setting.



FACTORY RESET WITH THE CONTROLLER ON

Press & Hold  for 10 seconds.

The controller will shut down and then turn back ON automatically.



INFORMATION

Once flashing temperature has stopped, the indicator displays the actual temperature measured at the core of the Mat.

Press and hold  or  to reach desired numbers faster.

Controller instructions can be found at any time at healthywavemat.com.

CONTROLLER INSTRUCTIONS

HEAT + PEMF CONTROLLER

POWER




Press the Power button to turn ON the controller. TEMP, TIMER, and PEMF indicators will light up on the screen.

TEMPERATURE

SET & ADJUST

Press the TEMP button once to turn ON the heat. Press TEMP again to initiate temperature adjustment mode (the temperature indicator will start flashing). Set the desired temperature and press TEMP again to save the setting.



°F or °C

UNITS SET

Press & Hold for 6 sec




To switch between Fahrenheit & Celsius.

TIMER

SET & ADJUST

Press the TIME button to initiate time adjustment mode (the time indicator will start flashing). Set the desired time and press the TIME button again to save the setting.



PEMF

SET & ADJUST

Press the PEMF button once to turn PEMF ON. Press PEMF again to initiate PEMF frequency adjustment mode (PEMF indicator will start flashing). Set the desired frequency (1-30 Hz) and press the PEMF button again to save the setting.



Press & Hold for 6 sec



To turn PEMF OFF.

FACTORY RESET

WITH THE CONTROLLER ON

Press & Hold  for 10 seconds.

The controller will shut down and then turn back ON automatically.



INFORMATION

Once flashing temperature has stopped, the indicator displays the actual temperature measured at the core of the Mat.

Press and hold  or  to reach desired numbers faster.

Default PEMF frequency is set to 7.83 Hz and displayed as 8 Hz.

PEMF can be turned on or off at any time with or without heat.

Controller instructions can be found at any time at healthywavemat.com.

CONTROLLER INSTRUCTIONS

HEAT + PHOTON CONTROLLER

POWER




Press the Power button to turn ON the controller. TEMP, TIMER, and PHOTON indicators will light up on the screen.

TEMPERATURE

SET & ADJUST

Press the TEMP button once to turn ON the heat. Press TEMP again to initiate temperature adjustment mode (the temperature indicator will start flashing). Set the desired temperature and press TEMP again to save the setting.



°F or °C

UNITS SET

Press & Hold
for 6 sec



To switch between
Fahrenheit & Celsius.

TIMER

SET & ADJUST

Press the TIME button to initiate time adjustment mode (the time indicator will start flashing). Set the desired time and press the TIME button again to save the setting.



PHOTON LIGHT

ON/OFF

Press the  button to turn lights ON or OFF.

FACTORY RESET

WITH THE CONTROLLER ON

Press & Hold  for 10 seconds.

The controller will shut down and then turn back ON automatically.



INFORMATION

Once flashing temperature has stopped, the indicator displays the actual temperature measured at the core of the Mat.

Press and hold  or  to reach desired numbers faster.

Photon can be turned on or off at any time with or without heat.

Controller instructions can be found at any time at healthywavemat.com.

CONTROLLER INSTRUCTIONS

HEAT + PEMF + PHOTON CONTROLLER


POWER

 Press the Power button to turn ON the controller. TEMP, TIMER, PEMF, and PHOTON indicators will light up on the screen.

TEMPERATURE

SET & ADJUST

Press the TEMP button once to turn ON the heat. Press TEMP again to initiate temperature adjustment mode (the temperature indicator will start flashing). Set the desired temperature and press TEMP again to save the setting.



°F or °C

UNITS SET

Press & Hold for 6 sec  To switch between Fahrenheit & Celsius.

TIMER

SET & ADJUST

Press the TIME button to initiate time adjustment mode (the time indicator will start flashing). Set the desired time and press the TIME button again to save the setting.



PEMF

SET & ADJUST

Press the PEMF button once to turn PEMF ON. Press PEMF again to initiate PEMF frequency adjustment mode (PEMF indicator will start flashing). Set the desired frequency (1-30 Hz) and press the PEMF button again to save the setting.



Press & Hold for 6 sec  To turn PEMF OFF.

PHOTON LIGHT

ON/OFF

Press the  button to turn lights ON or OFF.

FACTORY RESET

WITH THE CONTROLLER ON

Press & Hold  for 10 sec.

The controller will shut down and then turn back ON automatically.



INFORMATION

Once flashing temperature has stopped, the indicator displays the actual temperature measured at the core of the Mat.

Press and hold  or  to reach desired numbers faster.

Default PEMF frequency is set to 7.83 Hz and displayed as 8 Hz.

PEMF can be turned on or off at any time with or without heat.

For Rainbow and Chakra Series press "Photon Button" twice to light all seven colors.

Controller instructions can be found at any time at healthywavemat.com.

COMFORT

Each Healthy Wave mat is designed with the sole purpose of improving your health and wellness. It is a powerful wellness device. Every mat consists of several pounds of real gemstones. Our Mesh models utilize large round pieces of tourmaline and jade. These mats have a hard surface. Jet, TAO, TAJ, Rainbow Chakra, SOFT, Platinum, 360-Wrap Set,


as well as Portable models, utilize crushed gemstones that feel more like large pieces of sand. This is not the same as lying on a mattress. No matter what model you get, there will be varying degrees of firmness. This can be softened by covering your mat with additional layers of padding to cushion the surface.



JET CONTROLLER INSTRUCTIONS

HEAT + PEMF + PHOTON CONTROLLER

POWER



Press the I/O button to power up the controller (second press will power the controller down).




Press the Power button to turn the controller ON (the blue indicator will light up). Indicators TEMP, TIMER, PEMF, PHOTON will light up on the screen (second press will turn the controller OFF).

TEMPERATURE

SET & ADJUST

Press the TEMP button once to activate the heat (timer will be activated simultaneously). Press the TEMP button again to enter temperature adjustment mode (the temperature indicator will start flashing). Set the desired temperature and press the TEMP button again to save the setting.




Once the temperature indicator stops flashing, it will display the actual temperature measured at the core of the mat after 5 seconds.

TIMER

SET & ADJUST


Press the TIME button to enter timer adjustment mode (the time indicator will start flashing). Set the desired time (0:01, 0:02, ... 0:59, 1:00, 2:00, ... 12:00) and press the TIME button again to save the setting. The timer and heat will automatically turn off when the countdown reaches 0:00.



PEMF

SET & ADJUST

Press the PEMF button once to activate PEMF. 2nd press activates PEMF frequency adjustment mode (PEMF indicator will start flashing). Set the desired frequency (0.25, 0.5, 1, ... 30 Hz) and press the PEMF button again to save the setting. The 3rd press (after the adjustment mode exit due to any button press (including the PEMF button)) will deactivate the PEMF function.




PEMF is set as cycles of 20 min on and 100 min off until PEMF is turned OFF or TIMER countdown reaches 0:00.

PHOTON LIGHT

ON/OFF

Press the PHOTON button to turn the lights ON or OFF (the screen will display the countdown of the light operating time). Press & hold the TIME button to enter the lights operating time adjustment mode (the PHOTON indicator will start flashing). Set the time (1 to 9 hrs), and press the TIME button again to save the setting. The light will automatically turn off when the countdown reaches 0:00.



*Pressing any key will save the setting and exit adjustment mode; however, it is recommended to develop a habit of using buttons relevant to the function being adjusted.



CONTINUED

MEMORY MODE M1 - M4

Press the M1 to M4 buttons to activate a memorized mode. The screen will display the selected mode along with the corresponding temperature, time, PEMF, and photon settings (see factory presets below).

To adjust settings, follow the instructions provided above. To save new settings press & hold the corresponding M-button for 5 sec (a flashing blue indicator will confirm the successful update).



...



°F or °C UNITS SET

Press & hold the TEMP button for 6 sec to switch between Fahrenheit & Celsius.



FACTORY RESET WITH THE CONTROLLER ON

When the controller is turned ON, press & hold the Power button until the blue indicator light turns off. Release the Power button; the indicator light will then blink twice to signal that the controller has reset.



FACTORY SETTINGS

	Mode	Temp	Photon	PEMF	Timer
	Default	86°F (30°C)	1:00	7,83 Hz*	12:00
M1	Sleep	100°F (38°C)	off	3 Hz	12:00
M2	Relax	110°F (43°C)	off	8 Hz	2:00
M3	Energize	130°F (54°C)	0:30	12 Hz	1:00
M4	Recovery	140°F (60°C)	1:00	20 Hz	1:00

E1: Sensor is unplugged or connection lost.

E2: Sensor short circuit or mat overheating.

E3: Controller overheating.

Abn: PEMF coil or heating wire is disconnected. The indicator will flash for 5 sec, then vanish; all other functions will operate normally.

INFORMATION

- Once the temperature stops flashing, the indicator displays the actual temperature measured at the core of the mat.
- Upon activation, each function will operate based on the factory preset or the most recently saved settings.
- If there is no action for 5 seconds while in a function adjustment mode, the last displayed setting will be saved, and the adjustment mode will exit.
- Press and hold + or - to reach desired numbers faster.
- PEMF/PHOTON can be turned on or off at any time with or without heat.
- If there is no action for 5 min, the controller will dim and enter sleep mode (activated functions will continue to operate). Press any button to wake the screen.
- The controller is universal for 110V and 220V.

*Displayed as 8 Hz.

PLATINUM CONTROLLER INSTRUCTIONS

ADVANCED CONTROLLER MULTI-WAVE PEMF MAT



Power ON/OFF to activate controller.




Turns Infrared Heat ON or OFF. Indicator on screen confirms Infrared Heating is ON.




Indicator on screen "Hours ON" shows how much time is left for Infrared Heating. After this time the Infrared Heating will turn OFF automatically. "Actual Temperature" of the mat is shown on the screen.



To change "Set point Temperature" & "Hours On."



Press Heat Mode Button once to set Temperature. Set point temperature on screen will flash. Use UP/DOWN "H" buttons to increase or decrease the desired temperature set point.



Press Heat Mode Button again to set Timer for Infrared Heating. On Time indicator on screen will flash. Use UP/DOWN "H" buttons to increase or decrease the desired timer set point. (1 to 12 hours available).

Press Heat Mode Button again to save both settings. Both these new settings will appear on the screen.



Press to activate infrared sleep mode cycle. The Infrared Heating will be turned ON; Temperature set point will be 86°F (30°C), Heating ON time set to 12 hours.



Press to turn the photon lights ON or OFF. Red Lights will stay ON for 1 hour when turned ON and automatically turn OFF after 1 hour. Indicator on screen confirms lights are activated.



Turns standard PEMF cycle ON [10/15/20 min selectable] plus 100 min OFF cycle. (Do not press this when using any of the F1 to F12 PEMF programs.)


- When controller is first powered ON the default settings (as shown on screen) will be 20 minute ON cycle, "Sine" Wave Type, M (Medium) pulse duration, 50% Intensity.
- This standard PEMF cycle will automatically operate for 20 minutes ON time and 100 minutes OFF time. This will repeat for up to 12 hours if controller buttons are not touched after this. After 12 hours the controller will automatically turn OFF.
- During the PEMF 20 minute ON time the screen indicates "Pulse". During the 100 minute OFF time the screen indicates "Sleep".
- To stop the standard PEMF cycle at any time press the PEMF ON/OFF button.




To make changes to the standard PEMF cycle and save your settings.



Press once to make changes to PEMF ON time, "MIN/HR" will flash on the screen. Use UP/DOWN "P" buttons to increase or decrease the desired timer set point. (10, 15 or 20 minutes).



Press PEMF Mode Button again, "Wave Type will flash. Use UP/DOWN "P" buttons to increase or decrease the desired wave type (sine or square wave).



Press PEMF Mode Button again, "Frequency" will flash. Use UP/DOWN "P" buttons to increase or decrease the desired frequency (1 to 30 pulses per second [Hz]).



Press PEMF Mode Button again, "Pulse Duration" will flash. Use UP/DOWN "P" buttons to increase or decrease the desired PEMF Duration ((M) Medium, (F) Fast, (UF) Ultra Fast).



Press PEMF Mode Button again, "Intensity" will flash. Use UP/DOWN "P" buttons to increase or decrease the desired PEMF Intensity (10%, 25%, 50%, 75% , 100%).



CONTINUED



Button F1 to F10 turns ON/OFF PEMF Programs. Each PEMF Program is comprised of:

- 4 segments of 5 minute PEMF programmed selections. Total program runs for 20 Minutes PEMF ON time and 100 minutes OFF time. This will repeat for up to 12 hours if controller buttons are not touched after this. After 12 hours the controller will automatically turn OFF.
- Each PEMF 5 minute segment is programmed with a unique set of Wave Type/Frequency/Pulse Duration/Intensity.

TABLE 1 - PRESET PEMF FUNCTIONS DESCRIPTION

Program Button	Program Type	Default Frequency Selection (pulses per second)
F1	Low Frequencies	1, 3, 4, 6
F2	Medium Low Frequencies	7, 8, 10, 12
F3	Medium Frequencies	14, 15, 17, 18
F4	High Frequencies	19, 21, 23, 25
F5	Before Sleep	5, 4, 3, 2
F6	Pain Assist	15, 16, 19, 20
F7	Sport Injury & Strain Assist	24, 24, 25, 25
F8	General Regeneration	7.83, 7.83, 10, 10
F9	Earth Frequencies	7.83, 14, 21, 25
F10	Wake Up/Energy	7.83, 10, 11, 12

F11 F12 Button F1 to F10 turns ON/OFF PEMF Programs. Each PEMF Program is comprised of:



F11 F11 "Energize" 1 hour program: The following 3 - 20 minute programs will run in sequence, F10-F8-F6. The PEMF program ends after 1 hour and automatically turns OFF.




F12 F12 "Relaxation" 1 hour program: The following 3 - 20 minute programs will run in sequence, F9-F8-F5. The PEMF program ends after 1 hour and automatically turns OFF.

SCREEN INDICATORS DURING F1 TO F10 PROGRAMS:

- The selected program button has back light activated.
- Arrow flashes between A-B, B-C, C-D, D-A to indicate which 5 minute PEMF segment is active.
- "Pulse" shown during PEMF active time, "Sleep" during 100 Minute OFF time.
- Timer counts down for 5 minutes for each 5 minute PEMF segment.
- Wave Type/Frequency/Pulse Duration/Intensity is shown for each 5 minute PEMF segment.

HOW TO CREATE CUSTOM PROGRAMS



INFORMATION

Always plug this controller into a 120Volt surge protected power bar with a power ON/OFF switch. When the controller is not in use, turn the power bar switch OFF. The controller is for standard 120Volt operation only. (220 Volt available and will be shown on the controller label.)

For full instructions including custom programming please visit:
<https://healthywavemat.com/user-manuals-guides-ezp-6>





CARE

1. Protect your mat from moisture at all times. Cover the entire mat with a blanket, towel, or included waterproof cover (hereinafter - cover) to protect your body from direct heat and to protect the mat from sweat. You can also purchase a padded waterproof cover from healthywavemat.com.
2. For sanitary purposes change and wash covers after each user. You can also use a disposable cover for each use.
3. After use, set the mat to a high temperature to allow any moisture to evaporate. After 15-20 minutes unplug and power down the controller. Do not leave your mat unattended while it is powered on.
4. Clean the surface of the mat and controller with a dry cloth as needed or use a slightly damp sanitizing wipe.
5. Do not use chemicals or solvents during cleaning.
6. If you need to store your mat, use the dustproof bag provided.
7. To fold the mat, unplug the controller, and then fold the mat inward at predetermined places.
8. Keep the mat in a dry and ventilated place.
9. Your mat will last longer the less you fold or bend it.





NOTES

1. Our mats are good for people of all ages, including children. Pets can benefit from the mats, too!
2. We only sell new mats. The box you receive your package in may or may not contain our branding. The box may show signs of wear due to shipping.
3. Check our website for other products and accessories. We have special models designed for chairs, knees, arms, shoulders, and non-powered energy bedding.
4. Spread the word! Share your experience with friends, family, neighbors, social media, and groups. You may also become an affiliate or referrer.
5. For helpful videos, please visit [healthywavemat.com \[Support/Help\]](https://healthywavemat.com/Support/Help).



FEATURED THERAPIES

SET OF FEATURES DEPENDS ON MODELS

Our gemstone heating pads intended for muscular relaxation, temporary relief of stiffness, minor muscle and joint pain; for the temporary relief of joint pain associated with arthritis; for the temporary relief of muscle spasms, minor sprains and strains; for the temporary relief of minor muscular back pain; for the temporary relief of stiffness, and increase of circulation where applied.



HOT STONES

Healthy Wave Mat uses a variety of natural gemstones to create a natural healing experience. These gemstones release far-infrared rays and negative ions when heated. They can also impact pressure points to release tension in the body. Not only can they help improve physical discomfort, but also achieve the common goal of improved health.



FAR INFRARED

Advanced heating systems and natural gemstones deliver far infrared rays (FIR) deep into the body. FIR penetrates deep into the body, promoting the temporary relief of stiffness, joint pain, and elevating muscle relaxation. The higher you set the temperature, the greater the level of far-infrared rays. Even organic heat from your body activates the gemstones to deliver far-infrared benefits, but at a lower level.



NEGATIVE IONS

Negative ions are found in abundance in the mountains, near waterfalls and beaches. They are also naturally emitted from gemstones. When your mat is turned on, the negative ion output significantly increases. The higher you set the temperature, the greater the level of negative ions you receive. The output of negative ions varies from 300 to 3000 per cubic centimeter, depending on the type of gemstones and temperature.




PHOTON LIGHTS

The photon lights in Healthy Wave products use visible red light at a wavelength of 660 nm, that reaches just below the surface of your skin, penetrating the skin on a cellular level. For maximum benefits, photon lights should come into contact with the skin directly. Benefits are lessened when applied through a cloth or cover. Photon light therapy may help with skin rejuvenation and improvement of your overall levels of wellness.



PEMF

Our PEMF system simulates the Earth's electromagnetic field, which is delivered directly to your body. PEMF therapy is the process of supplying the body with much-needed magnetic fields, bypassing the electrical pollution that affects our bodies on a daily basis. By default, PEMF is emitted at the frequency of 7.83 Hz in a pulsed sinusoidal waveform - the most consistent with the biorhythms of the human body. The unique combination of revitalizing natural elements in our mats takes wellness to the next level.



UP TO **21**  LAYERS





OUR POLICIES

WARRANTY

Healthy Wave Mat offers a comprehensive manufacturers 15-Month Limited Warranty for all AC powered mats, and a 5-Year Limited Warranty for TAO, TAJ, Rainbow Chakra, SOFT, Mesh, Platinum, 360-Wrap Set, and any new series sold after August 1, 2020.

In order to make a warranty claim, you will need to provide a purchase confirmation that verifies the purchase condition like new, purchase date, model, and serial number, as well as the price paid (see page 17). For detailed information please visit <https://healthywavemat.com/shippinginfo>.

90-DAY TRIAL PERIOD

We are confident that you will love our products. That is why we offer a unique 45-Day Trial Period for any gemstone mat. You have 45 days to test, try, and experience the wonderful benefits of our unique product. If for any reason, you aren't satisfied within 45 days, you can contact our customer service and send it back to us; though, we fully expect you to love our product the day it arrives.

LIFETIME TRADE-IN

Healthy Wave Mat offers the world's best manufacturers lifetime trade-in policy. After the warranty period is over, regardless of your mat's condition, you can always send it to the manufacturer as a trade-in and use half of its original value for a new purchase. So, if something happens to your mat — it will never become worthless.



LIFETIME UPGRADE

If at any point you decide you want a larger or more advanced mat, you have the option to upgrade your current one. Please note that the new purchase must be at least twice the value of the mat you are intending to upgrade. As long as your mat is in good condition, simply send it in with assistance from our customer service team and you will get a credit to its full original value that will be put towards your new purchase of one or more items.





TROUBLESHOOTING

If for any reason your mat does not work, please be advised that this may be a minor issue that can be easily resolved. Make sure you set your mat correctly according to the instructions provided.

If there is an issue, try disconnecting the controller from the mat and the outlet; fold and unfold the mat to make sure that it lies flat so that it can function normally; turn on the controller again. Try using different sockets. Make sure you plug the controller all the way in on both ends and set the timer before trying to adjust the temperature (your mat will not heat up if the timer has not been set). Heat the covered mat at the maximum temperature for 30 minutes to check the controller is showing an increase in temperature. If you're still having problems, please prepare the following information and then contact us:

1. Have your order number or a copy of your invoice ready.
2. Determine if your product is still within the 1 to-5 year limited warranty period.
3. Find the model number and the serial number of your mat. They can be found on the tag of your mat.
4. Determine which controller you have (see pages 7 - 10) or visit www.healthywavemat.com [Support/Help section].
5. Your shipping address, including zip code.
6. Once you have this information ready, contact us by email at support@healthywavemat.com or call us at 360-233-2088.



DISCLAIMER

This information is not intended to cure, diagnose, advise, or treat any medical condition, nor is reading this a substitute for reading the product User Guide. Please consult with a physician before use of this product. These statements have not been evaluated by the Food and Drug Administration. These products are not intended to diagnose, treat, cure or prevent any disease. Information with regard to our products does not reflect or represent product claims. Healthy Wave Mat assumes no responsibility for the improper use of its products or for a user's self-diagnosis and/or treatment using its products. HealthyWaveMat products should not be used as a substitute for medically supervised therapy. If you suspect you suffer from clinical deficiencies, consult a licensed, qualified medical doctor. No action should be taken solely on the content provided by HealthyWavemat; regardless of any perceived scientific merit, readers should consult healthcare professionals on any matter related to their health. The information in HealthyWaveMat's materials was obtained from materials that are believed to be accurate, as presented by their respective authors, but HealthyWaveMat assumes no liability for any errors or personal interpretation. Users and customers who fail to consult their physicians prior to the purchase or use of any product assume the risk of any adverse effects incurred. The Food and Drug Administration has not approved or evaluated the information contained in the literature here or on the website. Any printed reference or direct link to a specific product is for your information or convenience only and may not be construed as an enticement to purchase and further is not intended or implied to be used in the mitigation, diagnosis, treatment, cure, or prevention of any disease. Except as specifically stated, neither HealthyWaveMat nor any of its directors, employees, or other representatives will be liable for damages arising out of or in connection with the use of this product or information. This is a comprehensive limitation of liability that applies to all damages of any kind, including (without limitation) compensatory, direct, indirect, or consequential damages, loss of data, income or profit, loss of or damage to property, and claims of third parties.

- All information in this booklet is subject to change without notice.
- HealthyWaveMat is not responsible for any mistakes or misprints.
- Check the latest updates at healthywavemat.com.



FAQ

Can I use it on the floor or bed?

Yes, our mats can be used on any flat surface such as a bed, massage table, or floor.

What size should I get to sleep on?

72" x 24" is recommended, but if you roll over during the night, you might consider getting 76" x 32" or 80" x 40". We also have Twin, Queen, or King sizes in some models.

How can I burn calories?

Using our mat at high settings gives your body similar effects to doing cardiovascular exercise. As a reference, an hour of cardio can burn up to 500 calories. Setting the infrared heat to a temperature of 122 F-140 F (50 C-60 C) will increase your heart rate, boost metabolism, increase your energy, and elevate blood circulation on a cellular level. You can also use the sauna effect for the calorie burner. Set the temperature to 158 F and cover your body with a regular or a mylar blanket. You will receive all the benefits of a sauna, including sweat and detoxification. At these temperatures, please keep it reasonably under one hour. And always drink plenty of water, before and after.

Is your mat protected against EMFs?

Yes, we use industry leading double insulated wiring for advanced EMF shielding. The system is then placed under EMF-filtering copper mesh layer. Mats are also made with grounding for additional EMF protection.

How long can I use the mat?

You can use this mat for short or extended periods – even overnight. It all depends on your preference. For extended or overnight use, settings should be at a low or comfortable setting.

If I have a serious illness, can I use this mat?

If you have any preexisting conditions that may be a cause for concern, or if you have heart disease, hypertension, or any other serious medical condition you should consult with a doctor before using this or any other medical device.

Can I use my mat overseas?

All mats use 110V power for the USA and Canada, but you can still use the mat overseas with a good power converter. Most models are now available in 220V power – just ask or look carefully through the website.

How is 20 lbs. of gemstones so affordable?

It would cost a fortune to buy small quantities. If you were to buy 100 lbs, it costs considerably less. At HealthyWaveMat, we buy thousands of pounds of certified gemstones each year. Because of this, we can afford to sell our mats at a fraction of its cost!

Can I wash and protect my mat?

Covering your mat fully with a towel, bed sheet, or waterproof cover will help avoid messy issues. If you must clean your mat, unplug it, and use a damp towel. If your mat is truly dirty, unplug it, and use a damp cloth with a small amount of soap. To rinse the soap, use a damp cloth and then let it air dry.

If I have more questions what is the best way to get answers?

You can always email us at support@healthywavemat.com, and we'll get right back to you, typically that same day.

What to do in case of a spill of something smelly?

1. Unplug and turn off the controller;
2. Use paper towels or dry soft cloth to absorb excess moisture as much as possible;
3. Thoroughly wipe the affected area with multi-surface cleaning wipes.
4. Spray the odor remover on the damaged area from 6-8 inches away. Leave the mat in a well-ventilated area until it dries completely.
5. Turn on the mat and set the maximum temperature. After the maximum is reached keep the mat running for 30-60 minutes. Do not cover the mat so that it can ventilate. Do not leave the mat unattended.
6. If it does not help, carefully rinse the affected area of the unconnected mat under the shower stream, minimizing the ingress of water on the adjacent areas.
7. Repeat steps 2-5 and blow-dry the mat from a reasonable distance until it is completely dry.



I have a slight headache after using the mat. What do I do?

A headache is a sign of detoxing too fast and being dehydrated. Use a lower temperature and drink more water. Enjoy the process. Maximize your wellness at a comfortable pace.

I see your mats are made with different gemstones. Which gemstone mat should I buy?

Each mat is made with a different combination of gemstones. What you choose should depend on ailments or the area of wellness you would like to focus on. In general, the best mats come with a combination of gemstones, along with built-in PEMF.

Can I fold the mat?

Yes. Although some models, such as ones with PEMF, are ridged and must be folded in predetermined areas. If you plan on storing your mat in a tight area, consider getting the Reverse or SOFT model. Keep in mind, the more you fold your mat, the faster it will lose durability.

HEALTHYWAVEMAT.COM

2026-2027



SCAN QR CODE TO SHOP

ONLINE

FAR INFRARED HEAT THERAPY MATS
FILLED WITH NATURAL GEMSTONES



925 Shepherd Avenue
New York, NY 10018



1-360-233-2088 (USA)
support@healthywavemat.com



@healthywavemats